



Representative photo(s): actual product specifications may vary from that shown in the photo

Positive Placement Spring-Wire Guide

AV SWG COMPONENT: .025" DIA

SKU / Article #: RW-04025

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 25

| | |
|---|---|
| Positive Placement Spring-Wire Guide SKU / Article #: RW-04025 | |
| Quantity | Component Description |
| 1 | Spring-Wire Guide/Tube Assembly: .025" (0.64 mm) dia. x 5-3/8" (14 cm) (Straight Soft Distal End) |

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

