

Femoral Artery Catheterization

ARTERIAL CATH SET: 18GA X 6 IN

SKU / Article #: SA-04018

Considerations



Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Femoral Artery Catheterization SKU / Article #: SA-04018	
Quantity	Component Description
1	Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
1	Scalpel: #11
1	Suture: 2-0 Silk with Curved Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

