

NO IMAGE  
AVAILABLE

Two-Lumen CVC

CVC SET: 2-LUMEN 4FR X 8CM

SKU / Article #: SA-15402

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.



Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

|  |  |
|--|--|
| Two-Lumen CVC<br>SKU / Article #: SA-15402 |  |
| Quantity                                   | Component Description  |
| 1  | Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 8 cm   |
| 1  | Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer |
| 1  | Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle  |
| 1  | Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW  |
| 1  | Syringe: 5 mL Luer-Slip  |
| 1  | Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm  |
| 1  | Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm  |
| 1  | SecondSite™ Adjustable Hub: Catheter Clamp   |
| 1  | SecondSite™ Adjustable Hub: Fastener   |
| 1  | Safety Scalpel: #11  |
| 2  | Valve: MicroCLAVE™ <sup>1</sup> Neutral Displacement Connector   |

<sup>1</sup> A trademark of ICU Medical, Inc.

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.