

NO IMAGE
AVAILABLE

Three-Lumen CVC

CVC SET: 3-LUMEN 5.5FR X 8CM

SKU / Article #: SA-15553

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.



Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC SKU / Article #: SA-15553	
Quantity	Component Description
1	Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 8 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Safety Scalpel: #11
3	Valve: MicroCLAVE™ ¹ Neutral Displacement Connector

¹ A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.