

## Three-Lumen CVC

CVC SET: 3-LUMEN 7FR X 20CM

SKU / Article #: SA-15703

### Highlights

---

Sharps Safety

### Considerations

---

Not made with natural rubber latex.

**Sales Quantity/Case: 10**

---

NO IMAGE  
AVAILABLE



*Representative photo(s): actual product specifications may vary from that shown in the photo*

Three-Lumen CVC  
SKU / Article #: SA-15703

Quantity	Component Description
1	Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
2	Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
1	Safety Scalpel: #11
2	Prep Sponge Swab
3	Valve: MicroCLAVE™ <sup>1</sup> Neutral Displacement Connector

<sup>1</sup> A trademark of ICU Medical, Inc.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

