



Three-Lumen CVC

CVC SET: 3-LUMEN 5.5FR X 13CM

SKU / Article #: SA-16553

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Three-Lumen CVC	
SKU / Article #: SA-16553	
Quantity	Component Description
1	Three-Lumen Catheter: 5.5 Fr. (1.9 mm OD) x 13 cm
1	Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Syringe: 5 mL Luer-Slip
1	Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
1	Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Safety Scalpel: #11
3	Valve: MicroCLAVE™ ¹ Neutral Displacement Connector

¹ A trademark of ICU Medical, Inc.

