

Two-Lumen CVC

CVC SET: 2-LUMEN 7FR X 20CM

SKU / Article #: SA-17702

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

NO IMAGE
AVAILABLE


Representative photo(s): actual product specifications may vary from that shown in the photo

Two-Lumen CVC
SKU / Article #: SA-17702

Quantity	Component Description
1	Two-Lumen Catheter: 7 Fr. (2.5 mm OD) x 20 cm
1	Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
1	Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
1	Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
1	SecondSite™ Adjustable Hub: Fastener
1	SecondSite™ Adjustable Hub: Catheter Clamp
1	Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
2	Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
1	Safety Scalpel: #11
2	Prep Sponge Swab
2	Valve: MicroCLAVE™ ¹ Neutral Displacement Connector

¹ A trademark of ICU Medical, Inc.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

