



Arterial Catheterization Set

ARTERIAL CATH. SET: 24 GA X 2.5CM

SKU / Article #: SAC-00324-PBX

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Arterial Catheterization Set	
SKU / Article #: SAC-00324-PBX	
Quantity	Component Description
1	Catheter: 24 Ga. x 2.5 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
1	Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle

Representative photo(s): actual product specifications may vary from that shown in the photo

