

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set

ARTERIAL CATH SET: 24 GA X 2.5CM

SKU / Article #: SAC-00324

Highlights

CE Mark

Considerations

Not made with natural rubber latex.


Sales Quantity/Case: 10

Arterial Catheterization Set

SKU / Article #: SAC-00324

Quantity	Component Description
1	Catheter: 24 Ga. x 2.5 cm Radiopaque Polyurethane
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
1	Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.