

*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Arterial Catheterization

### ARTERIAL CATH KIT: 20GA X 5CM

SKU / Article #: SAC-00520-CA

### Highlights

Sharps Safety

### Considerations

Not made with natural rubber latex.


### Sales Quantity/Case: 5

<b>Arterial Catheterization</b> <b>SKU / Article #: SAC-00520-CA</b>	
Quantity	Component Description
1	Catheter: 20 Ga. x 5 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 1-1/2" (4 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Syringe: 3 mL Luer-Lock
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist Sheet
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® <sup>2</sup> 10 cm x 12 cm
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of 3M Company.

## Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.