



## Arterial Catheterization Set

### ARTERIAL CATH SET: 24 GA X 5CM

SKU / Article #: SAC-00524

#### Highlights

CE Mark

#### Considerations


Not made with natural rubber latex.

**Sales Quantity/Case: 10**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arterial Catheterization Set	
SKU / Article #: SAC-00524	
Quantity	Component Description
1	Catheter: 24 Ga. x 5 cm Radiopaque Polyurethane
1	Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 9-3/4" (25 cm) (Straight Soft Tip on Both Ends)
1	Catheter: 24 Ga. x 3/4" (1.9 cm) Radiopaque over 26 Ga. Introducer Needle

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.