



Arterial Catheterization Set

ARTERIAL CATH SET: 18 GA X 8CM

SKU / Article #: SAC-00818-PBX

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set	
SKU / Article #: SAC-00818-PBX	
Quantity	Component Description
1	Catheter: 18 Ga. x 8 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 18 Ga. x 2" (5 cm) XTW

