



Arterial Catheterization Set

ARTERIAL CATH SET: 20GA X 12CM

SKU / Article #: SAC-01220-PBX

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization Set	
SKU / Article #: SAC-01220-PBX	
Quantity	Component Description
1	Catheter: 20 Ga. x 12 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW

