

## Arterial Catheterization

**SAC KIT: 20 GA X 12CM**

SKU / Article #: SAC-01220-X1A

### Highlights

Sharps Safety

### Considerations

Not made with natural rubber latex.  
Includes Lidocaine

**Sales Quantity/Case: 5**

*Representative photo(s): actual product specifications may vary from that shown in the photo*

Arterial Catheterization	
SKU / Article #: SAC-01220-X1A	
Quantity	Component Description
1	Catheter: 20 Ga. x 12 cm Radiopaque PEBA
1	Spring-Wire Guide, Marked: .021" (0.53 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on Both Ends)
1	Introducer Needle: 20 Ga. x 2-3/4" (7 cm) XTW
1	Injection Needle: SafetyGlide™ <sup>1</sup> 25 Ga. x 1" (2.54 cm)
1	Injection Needle: SafetyGlide™ <sup>1</sup> 23 Ga. x 1-1/2" (3.81 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

## Arterial Catheterization

SKU / Article #: SAC-01220-X1A

Quantity	Component Description
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
2	3 mL Applicator 2% CHG and 70% IPA ChloroPrep <sup>®2</sup> One-Step Solution with Hi-Lite Orange <sup>™</sup> Tint
1	SharpsAway <sup>®</sup> Disposal Cup
1	SharpsAway <sup>®</sup> II Locking Disposal Cup
1	Maximal Barrier Drape <sup>™</sup> with 4" fenestration
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist Sheet
1	Medication Label: 1% Lidocaine
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm <sup>®3</sup> 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Dressing: BIOPATCH <sup>®4</sup>
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle
1	HemoHopper <sup>®</sup> Fluid Receptacle






<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.

<sup>3</sup> A registered trademark of 3M Company.

<sup>4</sup> A registered trademark of Johnson & Johnson Corporation.

# Documents



SCROLL TO VIEW MORE →

The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.