



Arterial Catheterization

SAC KIT: 18 GA X 16CM

SKU / Article #: SAC-01618-CA

Highlights

Sharps Safety

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

Representative photo(s): actual product specifications may vary from that shown in the photo

Arterial Catheterization


SKU / Article #: SAC-01618-CA

Quantity	Component Description
1	Catheter: 18 Ga. x 16 cm Radiopaque PEBA
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
1	Injection Needle: SafetyGlide™ ¹ 25 Ga. x 1" (2.54 cm)
1	Injection Needle: SafetyGlide™ ¹ 23 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	Syringe: 5 mL Luer-Slip
1	SharpsAway® Disposal Cup
1	SharpsAway® II Locking Disposal Cup
1	Drape: 41" x 55" with 4" fenestration
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Needle Holder
1	Safety Scalpel: #11
1	Checklist Sheet
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Dressing: Tegaderm® ² 10 cm x 12 cm
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Bouffant Cap
1	Surgical Apparel: Mask with Eye Shield
1	Suture: 3-0 Silk with Curved Needle

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of 3M Company.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.