



Representative photo(s): actual product specifications may vary from that shown in the photo

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters

PSI SET: 8.5 FR

SKU / Article #: SI-09875-ESB

Highlights

CE Mark


Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters SKU / Article #: SI-09875-ESB	
Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.