



*Representative photo(s): actual product specifications may vary from that shown in the photo*

## Arrowg+ard Blue® Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

**PSI SET: 8.5 FR**

SKU / Article #: SI-09880-AG

### Highlights

Antimicrobial: Arrowg+ard Blue\*  
CE Mark

### Considerations

\*Chlorhexidine  
Not made with natural rubber latex.

**Sales Quantity/Case: 10**

Arrowg+ard Blue® Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters  
SKU / Article #: SI-09880-AG

Quantity	Component Description
1	Arrowg+ard Blue® Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.