



Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

PSI SET: 8.5 FR

SKU / Article #: SI-09880-SB

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 10


Representative photo(s): actual product specifications may vary from that shown in the photo

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

SKU / Article #: SI-09880-SB

Quantity	Component Description
1	Arrow-Flex® Sheath: 8.5 Fr. (2.8 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.