



## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 - 8 Fr. Catheters

PSI SET: 9 FR

SKU / Article #: SI-09903-E

### Highlights

CE Mark

### Considerations

Not made with natural rubber latex.

**Sales Quantity/Case: 10**



*Representative photo(s): actual product specifications may vary from that shown in the photo*

### Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 - 8 Fr. Catheters

SKU / Article #: SI-09903-E

Quantity	Component Description
1	Arrow-Flex® Sheath: 9 Fr. (2.9 mm tip ID) x 10 cm and Tissue Dilator
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
1	Dust Cap: Non-Vented
1	Stopcock: 4-way High-Flow
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

# Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.