



Arrowg+ard Blue® MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 7.5 Fr. Catheters

MULTI-LUMEN/PSI SET

SKU / Article #: SI-21142

Highlights

Antimicrobial: Arrowg+ard Blue*
CE Mark

Considerations

*Chlorhexidine
Not made with natural rubber latex.

Sales Quantity/Case: 5


Representative photo(s): actual product specifications may vary from that shown in the photo

Arrowg+ard Blue® MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve
for use with 7 - 7.5 Fr. Catheters
SKU / Article #: SI-21142

Quantity	Component Description
1	Two-Lumen Central Venous Access Device: 4" (10 cm) Radiopaque Polyurethane with 9 Fr. Distal Lumen, Arrowg+ard® Antimicrobial Surface Treatment ¹ , Integral Hemostasis Valve, Extension Line Clamps
1	Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
1	Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
1	Pressure Transduction Probe
1	Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
1	Tissue Dilator: 9 Fr. (3.0 mm) x 20.7 cm
1	Obturator: 8 Fr.
1	Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
1	Dust Cap: Vented
3	Dust Cap: Non-Vented
1	Stopcock: 4-way High-Flow
3	Gauze Pad: 4" x 4" (10 cm x 10 cm)

¹ Licensed under US Patent No. 6,706,024.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

