

Access Kit

ACCESS KIT

SKU / Article #: UK-25000-MSB

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

Access Kit

SKU / Article #: UK-25000-MSB

Quantity	Component Description
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Injection Needle: Eclipse™ ¹ 23 Ga. x 1" (2.54 cm)
1	Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm)
1	Syringe: 5 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	Syringe: 10 mL Luer-Lock
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation
1	Towel
1	Needle Holder
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Insertion Checklist
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
6	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Cap
1	Dressing: OpSite® ² 4-3/4" x 4" (12 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Mask with Shield
1	Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)
1	Suture: (0) Silk with Curved Needle
1	Tray: Prep

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of T.J. Smith & Nephew, Limited.

