

Access Kit

ACCESS KIT

SKU / Article #: UK-25000-MSB

Highlights

CE Mark

Considerations

Not made with natural rubber latex.

Sales Quantity/Case: 5

NO IMAGE
AVAILABLE

Representative photo(s): actual product specifications may vary from that shown in the photo

Access Kit



SKU / Article #: UK-25000-MSB

Quantity	Component Description
1	Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
1	Injection Needle: Eclipse™ ¹ 23 Ga. x 1" (2.54 cm)
1	Injection Needle: Eclipse™ ¹ 25 Ga. x 1" (2.54 cm)
1	Syringe: 5 mL Luer-Lock
1	Syringe: 3 mL Luer-Lock
2	Syringe: 10 mL Luer-Lock
1	SharpsAway® II Locking Disposal Cup
1	SharpsAway® Disposal Cup
1	Towel: 13" x 18" (33 cm x 46 cm)
1	Maximal Barrier Drape: 55" x 94-1/2" (140 cm x 240 cm) with 4" (10 cm) fenestration, with adhesive and covered perforation
1	Needle Holder
1	Safety Scalpel: #11
1	Sterile Procedure Sign
1	Insertion Checklist
2	Gauze Pad: 2" x 2" (5 cm x 5 cm)
6	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Surgical Apparel: Cap
1	Dressing: OpSite® ² 4-3/4" x 4" (12 cm x 10 cm)
1	Surgical Apparel: Impervious Gown
1	Surgical Apparel: Mask with Shield
1	Transducer Cover: 4-12/16" x 95" (12 cm x 244 cm)
1	Suture: (0) Silk with Curved Needle
1	Tray: Prep

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of T.J. Smith & Nephew, Limited.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.