

Femoral Arterial Line Catheterization

ACCESS TRAY

SKU / Article #: UM-04018X

Considerations

Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 10

NO IMAGE
AVAILABLE


Representative photo(s): actual product specifications may vary from that shown in the photo

Femoral Arterial Line Catheterization

SKU / Article #: UM-04018X

Quantity	Component Description
1	Indwelling Catheter: 18 Ga. x 16 cm with Integral Suture Wing
1	Spring-Wire Guide: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
1	Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
1	Introducer Needle: 18 Ga. x 3" (7.62 cm) RW
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 10 mL Luer-Lock
2	Syringe: 5 mL Luer-Slip
2	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	Catheter Clamp
1	Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
1	Scalpel: #11
1	Prep Razor
2	Prep Sponge Swab
4	Gauze Pad: 2" x 2" (5 cm x 5 cm)
8	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle
1	Tray: Prep

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.