



Representative photo(s): actual product specifications may vary from that shown in the photo

Radial Artery Catheterization

RA CATH KIT: 20 GA X 1-3/4"

SKU / Article #: WBH-04020-1

Highlights

Sharps Safety

Considerations


Not made with natural rubber latex.
Includes Lidocaine

Sales Quantity/Case: 10

Radial Artery Catheterization SKU / Article #: WBH-04020-1	
Quantity	Component Description
1	Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
1	Injection Needle: 25 Ga. x 1" (2.54 cm)
1	Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
1	Syringe: 3 mL Luer-Lock
1	Skin Protectant Prep
1	5 mL 1% Lidocaine HCl Solution and Alcohol Prep
1	1.5 mL Applicator 2% CHG and 70% IPA ChloraPrep ^{®1} One-Step Frepp [®] Solution
1	Wing Clip: Suture
1	SharpsAway [®] II Locking Disposal Cup
1	Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
1	Safety Scalpel: #11
4	Gauze Pad: 4" x 4" (10 cm x 10 cm)
1	Suture: 3-0 Braided Silk with Straight Needle

¹ A registered trademark of CareFusion or one of its subsidiaries.

Documents



The products in this catalog may not be available in all countries. Please contact your local representative.

Indications, contraindications, warnings and instructions for use can be found in the Product Instructions for Use available under the Documents section on this page, where an IFU is available.

