

ARROW® EZ-IO®
INTRAOSSUEOUS VASCULAR ACCESS

Manual Needle Set Training Kit

38.5 mm 15 ga

**38.5 mm x 1.8 mm
40 kg or over**

Contents

Qty.

EZ-IO® 38.5 mm 15 ga MANUAL NEEDLE SET	6
EZ-Stabilizer® Dressing	6
EZ-Connect® Extension Set (Priming Volume – Approx. 1.0 mL)	6
Quick Reference Guide	6
Sharps Block	6
EZ-IO® Wristband	6
Instructions for Use	1

NO STERNUM



8262 Rev. 02

Teleflex®

 Manufactured for:
Teleflex Medical
IDA Business &
Technology Park,
Dublin Rd, Athlone,
Co. Westmeath, Ireland
Teleflex.com/military
Ph: 1.866.479.8500

REF
9094TK

U.S. and foreign patents issued and pending.

Product of the United States, Packaged in the United States.

NOT FOR HUMAN USE



Fluid path of needle and extension set is
non-sterile, non-DEHP, and not
made with natural rubber latex.