


38.5 mm 15 ga

**38.5 mm x 1.8 mm
40 kg or over**

Contents	Qty.
EZ-IO® 38.5 mm 15 ga MANUAL NEEDLE SET†	6
EZ-Stabilizer® Dressing	6
EZ-Connect® Extension Set (Priming Volume – Approx. 1.0 mL)	6
Quick Reference Guide	6
Sharps Block	6
EZ-IO® Wristband	6
Instructions for Use	1

NO STERNUM



8293A, Rev. 00 (2021-09)

Teleflex®

	
	
<p>U.S. and foreign patents issued and pending. Product of the United States, Packaged in the United States.</p>	
	
<p>Rx only. Fluid path of needle and extension set is sterile, non-pyrogenic, non-DEHP, and not made with natural rubber latex. †MR unsafe.</p>	