



Needle + Stabilizer Kit

REF 9079P



45 mm 15 ga

45 mm X 1.8 mm
40 kg or over

ARROW® EZ-IO®
INTRASOSSEOUS VASCULAR ACCESS

CONTENTS	QTY
EZ-IO® 45 mm 15 ga Intraosseous Needle Set	1
EZ-Stabilizer™ Dressing	1
EZ-Connect® Extension Set 1 (Priming Volume - Approx. 1.0 ml)	1
EZ-IO® Wristband	1
NeedleWISE® 1-Port Sharps Block	1

Fluid path of needle and extension set is sterile, non-pyrogenic, non-DEHP, and not made with natural rubber latex.

WARNING: This product contains a chemical known to the State of California to cause cancer, birth defects, or other reproductive harm.

CE
0086

STERILE EO



Rx ONLY



Manufactured for:
Teleflex Medical
IDA Business & Technology Park,
Dublin Rd, Athlone,
Co. Westmeath, Ireland
ArrowEZIO.com
Ph: 1.866.479.8500

Product of the U.S.A.

8772 Rev 02

Teleflex®

US and foreign patents
issued and pending.

Not For Sternal Use
PEEL HERE

PMS 281

PMS 423

PMS 107

PMS Black