



# ARROW<sup>+</sup> ARROWgard Blue<sup>®</sup> MAC<sup>™</sup> Two-Lumen Central Venous Access Kit Cath-Gard<sup>®</sup> for use with 7.5 - 8 Fr. Catheters



REF **CDC-21242-1A**

**9** Fr. distal lumen

**2** Lumen

**11.5** cm device length

**.035** inch dia. spring-wire guide

**ErgoPack<sup>®</sup>**

**Contents:**

- 1: Two-Lumen Central Venous Access Device<sup>1</sup>: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, ARROWgard<sup>®</sup> Antimicrobial Surface Treatment<sup>2</sup>, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide<sup>™</sup> 25 Ga. x 1" (2.54 cm)
- 1: Local Anesthetic Component
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide<sup>™</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock

- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep<sup>®4</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Cath-Gard<sup>®</sup> Catheter Contamination Shield: 80 cm with TwistLock<sup>™</sup> Adapter
- 1: Dust Cap: Vented
- 1: Dust Cap: Non-Vented
- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow

- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Surgical Apparel: Impervious Gown
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Tubing<sup>†</sup>: 20" (50 cm) Extension
- 1: HemoHopper<sup>®</sup> Fluid Receptacle

<sup>1</sup>Licensed under US Patent No. 5,250,038.  
<sup>2</sup>Licensed under US Patent No. 6,706,024.  
<sup>3</sup>A trademark of Becton, Dickinson and Company.  
<sup>4</sup>A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek<sup>®</sup> packaging



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine acetate, silver sulfadiazine and/or sulfa drugs.**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.


Not made with natural rubber latex

Precaution: Contains Phthalate: DEHP.

Store between 20 - 25°C (68 - 77°F).

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.70	29920
Distal (9 Fr.) with 8 Fr. catheter	1.20	12180
Distal (9 Fr.) with 7 Fr. catheter	1.55	19810
Proximal (12 Ga.) no catheter	0.55	9340

\* Priming volumes are approximate and are done without accessories.  
† Based on 100 cm head height.



Manufacturer:  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex<sup>®</sup>**  
[www.teleflex.com](http://www.teleflex.com)

C-21242-200E, Rev. 01 (2020-12)



# ARROW<sup>+</sup> ARROWgard Blue<sup>®</sup> MAC<sup>™</sup> Two-Lumen Central Venous Access Kit Cath-Gard<sup>®</sup> for use with 7.5 - 8 Fr. Catheters



REF **CDC-21242-1A**

**ErgoPack<sup>®</sup>**

**9** Fr. distal lumen

**2** Lumen

**11.5** cm device length

**.035** inch dia. spring-wire guide