



Arterial Catheterization Kit

REF NA-04550-K

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Contents:

- | | |
|--|---|
| <ul style="list-style-type: none"> 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque Polyurethane with Integral Suture Wing 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other) 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW 1: Injection Needle: SafetyGlide™¹ 25 Ga. x 1" (2.54 cm) 1: Injection Needle: SafetyGlide™¹ 22 Ga. x 1-1/2" (3.81 cm) RW 1: Syringe: 3 mL Luer-Lock | <ul style="list-style-type: none"> 2: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®² One-Step Solution with Hi-Lite Orange™ Tint 1: SharpsAway II™ Locking Disposal Cup 1: Drape: 41" x 55" with 4" fenestration 1: Towel 1: Scalpel: #11 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) 4: Gauze Pad: 4" x 4" (10 cm x 10 cm) <p>¹A trademark of Becton, Dickinson and Company.
 ²A registered trademark of CareFusion or one of its subsidiaries.</p> |
|--|---|



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex



Store between 15 - 30°C (59 - 86°F).
 Avoid freezing and excessive heat above 40°C (104°F)

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

DuPont™ Tyvek®
packaging

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®

www.teleflex.com

IA-04550-100B (2/17)



Arterial Catheterization Kit

REF NA-04550-K

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide