

REF CA-04706

# Central Venous Catheterization Set

# ARROW®

Catheter Diameter: 14 Ga., 2.2 mm

Priming Volume: 0.27 mL

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

Flow Rate: 7000 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

**Contents:**

- One: . . . . . Indwelling Catheter: 14 Ga. x 6" (16 cm) Radiopaque Polyurethane with Side Holes, Integral Suture Wing
- One: . . . . . Spring-Wire Guide: .032" (.81 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- One: . . . . . Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- One: . . . . . Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- One: . . . . . Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- One: . . . . . Tissue Dilator
- One: . . . . . T-Port Connector
- One: . . . . . Scalpel: #11

**Rx only.**


**Warning: Prior to use read all package insert warnings, precautions, and instructions.**

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide. Sharp objects included in this set must be discarded in a safe way.

This product contains dry natural rubber.



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA



1A-04706-1006 (11/09)



REF CA-04706 14 Ga. 16 cm catheter length Central Venous Catheterization Set

.032 inch dia. spring-wire guide

Peel to Open