

**One-Piece
Sheath
Design**

REF CA-09801

**Percutaneous Sheath
Introducer Kit
with Integral Hemostasis Valve/Side Port
for use with 7 – 7.5 Fr. Catheters**

ARROW®

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm)
Arrow-Flex® Radiopaque Polyurethane with Integral
Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia.
x 17-3/4" (45 cm) (Straight Soft Tip on One End –
"J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 1-9/16" (3.97 cm)
Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Pressure Transduction Probe
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Introducer Needle: 18 Ga. x 2-1/2"
(6.35 cm) XTW and 5 mL Arrow® Raulerson
Spring-Wire Introducer Syringe
- 1: Cath-Gard® Catheter Contamination Shield:
80 cm with Sterile Tape
- 1: Obturator: 8 Fr.
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: CSR Wrap
- 1: Drape: 24" x 36"
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Scalpel: #11
- 1: Stopcock: 3-way High-Flow
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 2: Prep Sponge Swabs
- 1: Suture: 3-0 Nylon with Straight Needle
- 1: Tray: Prep


Rx only.

**Warning: Read all package insert warnings, precautions,
and instructions prior to use. Failure to do so may result
in severe patient injury or death.**

Sterile. Single use — do not resterilize. Sterilized by
ethylene oxide. Do not use if package has been previously
opened or damaged. Fluid paths components are
non-pyrogenic. Contains no medication.

This Product Contains Dry Natural Rubber.

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA



IA-09801-100L (1/12)

ARROW®
INTERNATIONAL

REF CA-09801 8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit

Peel to
Open