

REF CA-09903-LF

LATEX-FREE

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7.5 - 8 Fr. Catheters

ARROW®

Contents:

- One: Sheath: 9 Fr. x 4" (10 cm)
Arrow-Flex® Radiopaque Polyurethane
with Integral Hemostasis Valve/Side Port
and Tissue Dilator
- One: Spring-Wire Guide, Marked: .035"
(.89 mm) dia. x 17-13/16" (45 cm)
(Straight Soft Tip on One End –
"J" Tip on Other)
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm)
Radiopaque over 20 Ga. RW
Introducer Needle
- One: Injection Needle: 25 Ga. x 1" (2.54 cm)
and 3 mL Luer-Lock Syringe
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
and 5 mL Luer-Slip Syringe
- One: Introducer Needle: 18 Ga. x 2-1/2"
(6.35 cm) XTW and 5 mL Luer-Slip Syringe
- One: Dust Cap: Non-Vented
- One: Cath-Gard® Catheter Contamination Shield:
80 cm with TwistLock™ Adapter
- One: Drape: 24" x 36" with 4" fenestration
- One: CSR Wrap
- Two: Gauze Pads: 2" x 2"
- Five: Gauze Pads: 4" x 4"
- One: Obturator: 8 Fr.
- One: Scalpel: #11
- One: Stopcock: 4-way High-Flow
- One: Suture: 3-0 Silk with Straight Needle
- One: SharpsAway® Disposal Cup

Rx only.


Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA



1A-09903-102D (4/09)

ARROW®

INTERNATIONAL

REF
CA-09903-LF

9 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit

Peel to Open
LATEX-FREE