



MAC™ Two-Lumen Central Venous Access Kit with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

REF CA-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Central Venous Access Device¹: 4" (10 cm) Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 20.7 cm
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
- 3: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Drape: 24" x 36" (61 cm x 91 cm)
- 1: Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

- fr** Kit d'accès veineux central à deux lumières MAC avec valve hémostatique intégrée, à utiliser avec les cathéters
- de** Einführbesteck für zweilumigen, MAC-zentralen Venenkatheter mit Hämostaseventil mit Kathetern der Größen 7 - 8 Fr.
- it** Kit di accesso venoso centrale a due lumi MAC con valvola emostatica incorporata per l'uso con cateteri da 7 - 8 Fr.
- pt** Kit para Acesso Venoso Central MAC com Lúmen Duplo, Válvula de Hemostase Integral, para utilizar com Cateters de 7 - 8 Fr.
- es** Juego para acceso venoso central de dos luces MAC con válvula hemostática integral para uso con catéteres de entre 7 y 8 Fr
- sv** MAC tvåkanaligt uppsättning för centralt ventillträde med integrerad hemostasventil för användning med 7 - 8 Fr katetrar

¹Licensed under US Patent No. 5,250,038.



Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.65	33000
Distal (9 Fr.) with 8 Fr. catheter	0.90	10500
Proximal (12 Ga.) no catheter	0.60	13000

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.

EC REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IA-11142-130C (02/20)



MAC™ Two-Lumen Central Venous Access Kit with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

REF CA-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide