

Multi-Lumen CVC Kit




REF ASK-15703-PCA

7^{Fr.} | 3^{Lumen} | 20^{cm} catheter length | .032^{inch dia.} spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: Eclipse™¹ 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
- 1: Syringe: 10 mL Luer-Lock
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 68" x 110" with 4" fenestration
- 1: Towel
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Flow Rate Information Card with Injection Log
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®² CHG 3-1/2" x 4-1/2" (8.5 cm x 11.5 cm)**
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper® Fluid Receptacle

¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of 3M Company.



| Lumen | Priming Volume* (mL) | Flow Rate† (mL/hr) | MAX Pressure Injection Flow Rate‡ (mL/sec) |
|-------------------|----------------------|--------------------|--|
| Distal (16 Ga.) | 0.44 | 3200 | 10 |
| Medial (18 Ga.) | 0.39 | 1500 | 5 |
| Proximal (18 Ga.) | 0.39 | 1600 | 5 |

* Priming volumes are approximate and are done without an injection cap.
 † Based on 100 cm head height.
 ‡ Pressure injection flow rates are determined using media of 11.8 centipoise viscosity, with 152 cm administration set tubing.

~~LATEX~~
 Not made with natural rubber latex