



# ARROW®

## PICC Peripherally Inserted Central Venous Catheter with Placement Wire

REF **PC-01351-TW**

**3** Fr. | **55** cm catheter length

**Contents:**

- 1: Indwelling Catheter, Trimmable: 3 Fr. x 55 cm Radiopaque Polyurethane with T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire
- 1: Dust Cap: Non-Vented
- 1: Fastener: Catheter Clamp
- 1: Catheter Trimmer
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure


- fr** Cathéter central veineux à insertion périphérique avec fil de placement
- de** Peripher eingeführter, zentraler Venenkatheter mit Platzierungsdraht
- it** Catetere per vena centrale a inserimento periferico con guida metallica di posizionamento
- pt** Cateter Venoso Central de Introdução Periférica com Fio de Colocação
- es** Catéter venoso central para inserción periférica con dispositivo de colocación metálico
- sv** Perifert införd central venkater med placeringsledare

All components are CE 0086 unless otherwise noted.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (20 Ga.)	0.30	90	600

\* Priming volumes are approximate and are done without accessories.  
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.  
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig.



**Manufacturer:**  
**Arrow International, Inc.**  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IC-01351-131C (11/17)



## PICC Peripherally Inserted Central Venous Catheter with Placement Wire

# ARROW®

REF **PC-01351-TW**

**3** Fr. | **55** cm catheter length