



PICC with Blue FlexTip®

ARROW®

REF **PC-01451**

4 Fr. | **55** cm
catheter
length

Contents:

- 1: Indwelling Catheter: 4 Fr. x 55 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp,
- Integral Suture Wing, Contamination Guard
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Fastener: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure

DuPont™ Tyvek®
insert

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65


Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: This device is contraindicated wherever there is presence of device related infections, previous or current thrombosis. Clinical assessment of patient must be completed to ensure no contraindications exists.



Not made with
natural rubber latex



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.50	850	3700

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

Manufacturer:

Arrow International LLC

Subsidiary of Teleflex Incorporated

3015 Carrington Mill Blvd.

Morrisville, NC 27560 USA

Teleflex®

www.teleflex.com

IC-01451-100H (2021-04)



PICC with Blue FlexTip®

REF **PC-01451**

4 Fr. | **55** cm
catheter
length

ARROW®