



# ARROW®

## PICC with Blue FlexTip® and Placement Wire

REF **PC-01451-W**

**4** Fr. | **55** cm catheter length

**Contents:**

- 1: Indwelling Catheter: 4 Fr. x 55 cm Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Fastener: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure

DuPont™ Tyvek®  
material

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).


Fluid path components are non-pyrogenic.



STERILE EO



Not made with natural rubber latex



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.50	850	3700

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

Catheter may be inserted through an Arrow® 14 Ga. Peel-Away introducer Sheath.

Manufacturer:

Arrow International, Inc.

Subsidiary of Teleflex Incorporated

2400 Bernville Road

Reading, PA 19605 USA

**Teleflex®**

[www.teleflex.com](http://www.teleflex.com)

IC-01451-101F (9/17)



## PICC with Blue FlexTip® and Placement Wire

REF **PC-01451-W**

**4** Fr. | **55** cm catheter length

# ARROW