



# Trimable PICC with Placement Wire

ARROW®

REF **PC-01451-TW**

**4** Fr. | **55** cm catheter length

**Contents:**

- 1: Indwelling Catheter, Trimmable: 4 Fr. x 55 cm Radiopaque Polyurethane with T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire
- 1: Dust Cap: Non-Vented
- 1: Fastener: Catheter Clamp
- 1: Catheter Trimmer
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure

DuPont™ Tyvek® | made in USA

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Contraindications: This device is contraindicated wherever there is presence of device related infections, previous or current thrombosis. Clinical assessment of patient must be completed to ensure no contraindications exists.**




STERILE EO



~~LATEX~~

Not made with natural rubber latex



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.50	850	3700

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

Above values apply to the untrimmed catheter. Trimming the catheter will alter priming volumes and flow rates.

**Manufacturer:**

Arrow International, Inc.

Subsidiary of Teleflex Incorporated

2400 Bernville Road

Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

IC-01451-102F (7/17)



# Trimable PICC with Placement Wire

ARROW®

REF **PC-01451-TW**

**4** Fr. | **55** cm catheter length