



Two-Lumen PICC with Placement Wire

ARROW®

REF **PC-01552-W**

5 Fr. | **2** Lumen | **55** cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 5 Fr. x 55 cm Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line
- Clamps, Contamination Guard and Placement Wire
- 2: Dust Cap: Non-Vented
- 1: Fastener: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure

DuPont™ Tyvek®
material

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Not made with natural rubber latex



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.40	575	2450
Proximal (20 Ga.)	0.40	160	1025

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

Manufacturer:

Arrow International, Inc.

Subsidiary of Teleflex Incorporated

2400 Bernville Road

Reading, PA 19605 USA

Teleflex®

www.teleflex.com

IC-01552-101E (9/17)



Two-Lumen PICC with Placement Wire

ARROW

REF **PC-01552-W**

5 Fr. | **2** Lumen | **55** cm catheter length