



Trimable PICC with Placement Wire

ARROW®

REF PC-05531-TW

3 Fr. | 55 cm catheter length

Contents:

- 1: Indwelling Catheter, Trimmable: 3 Fr. x 55 cm Radiopaque Polyurethane with T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire
- 1: Dust Cap: Non-Vented
- 1: Fastener: Catheter Clamp
- 1: Catheter Trimmer
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Paper Tape Measure

DuPont™ Tyvek®
material

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: This device is contraindicated wherever there is presence of device related infections, previous or current thrombosis. Clinical assessment of patient must be completed to ensure no contraindications exists.




STERILE EO



~~LATEX~~

Not made with natural rubber latex



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (20 Ga.)	0.30	90	600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

‡ Pump flow rates determined at maximum pump pressure of 10 psig.

Above values apply to the untrimmed catheter. Trimming the catheter will alter priming volumes and flow rates.

Manufacturer:

Arrow International, Inc.

Subsidiary of Teleflex Incorporated

2400 Bernville Road

Reading, PA 19605 USA

Teleflex®

www.teleflex.com

IC-05531-100C (7/17)



Trimable PICC with Placement Wire

ARROW®

REF PC-05531-TW

3 Fr. | 55 cm catheter length