

Product No. IC-09880

Percutaneous Sheath Introducer Set

with Integral Hemostasis Valve/Side Port
for use with 7 - 7.5 Fr. Catheters

ARROW®

Special Order Product

Contents:

- One: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex®
Radiopaque Polyurethane with Integral
Hemostasis Valve/Side Port and Tissue Dilator
- One: Spring-Wire Guide: .035" (.89 mm) dia.
x 17-13/16" (45 cm) (Straight Soft Tip
on One End – "J" Tip on Other)
- One: Introducer Needle: 18 Ga. x
2-1/2" (6.35 cm) XTW
- One: Cath-Gard® Catheter Contamination Shield:
80 cm with Sterile Tape
- One: Obturator: 8 Fr.
- One: Dust Cap: Non-Vented
- One: Safety Scalpel: #11
- One: Stopcock: 3-way High-Flow
- Three: Gauze Pads: 4" x 4"

Sheath O.D.: 3.6 mm

Rx only.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains dry natural rubber.



Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

IC-09880-103F, Rev. 1 (9/10)

ARROW
INTERNATIONAL
Special Order Product

Product No.
IC-09880

8.5 Fr.

10 cm
sheath
length

.035 inch dia.
spring-wire
guide

Percutaneous Sheath Introducer Set

Peel to
Open