



Two-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-12702**

7 Fr. | 2 Lumen | 16 cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.




- fr** Cathéter veineux central à deux lumières avec Blue FlexTip
- de** Zweilumiger, zentraler Venenkatheter mit Blue-FlexTip
- it** Catetere per vena centrale a due lumi con Blue FlexTip
- pt** Cateter Venoso Central com Lúmen Duplo e Blue FlexTip
- es** Catéter para vena central de dos luces con Blue FlexTip
- sv** Tvåkanalig centralvenös kateter med Blue FlexTip

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.54	5300
Proximal (18 Ga.)	0.37	2100

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Teleflex Medical
 IDA Business and Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



IC-12702-135A (03/19)



Two-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-12702**

7 Fr. | 2 Lumen | 16 cm catheter length