



## Arrow-Howes™ Multi-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-14703**

7 Fr. | 3 Lumen | 30 cm catheter length

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.



- fr** Cathéter veineux central à lumières multiples Arrow-Howes avec Blue FlexTip
- de** Viellumiger, zentraler Arrow-Howes-Venenkatheter mit Blue-FlexTip
- it** Catetere per vena centrale multilume Arrow-Howes con Blue FlexTip
- pt** Cateter Venoso Central Arrow-Howes com Lúmen Múltiplo e Blue FlexTip
- es** Catéter para vena central Arrow-Howes de luz múltiple con Blue FlexTip
- sv** Arrow-Howes flerkanalig centralvenös kateter med Blue FlexTip


Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm**  
[www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.54	2440
Medial (18 Ga.)	0.45	1080
Proximal (18 Ga.)	0.48	1150

\* Priming volumes are approximate and are done without accessories.  
† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



**Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**CE** EC REP 0086 Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IC-14703-135A (06/19)



## Arrow-Howes™ Multi-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-14703**

7 Fr. | 3 Lumen | 30 cm catheter length