



Arrow-Howes™ Multi-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-15703**

7 Fr. | 3 Lumen | 20 cm catheter length

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.




- fr** Cathéter veineux central à lumières multiples Arrow-Howes avec Blue FlexTip
- de** Viellumiger, zentraler Arrow-Howes-Venenkatheter mit Blue-FlexTip
- it** Catetere per vena centrale multilume Arrow-Howes con Blue FlexTip
- pt** Cateter Venoso Central Arrow-Howes com Lúmen Múltiplo e Blue FlexTip
- es** Catéter para vena central Arrow-Howes de luz múltiple con Blue FlexTip
- sv** Arrow-Howes flerkanalig centralvenös kateter med Blue FlexTip

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov

Fluid path components are non-pyrogenic.



Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

CE EC REP 0086
Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Teleflex®
www.teleflex.com

IC-15703-135A (03/19)



Arrow-Howes™ Multi-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-15703**

7 Fr. | 3 Lumen | 20 cm catheter length