



Two-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-17702**

7 Fr. | **2** Lumen | **20** cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.




- fr** Cathéter veineux central à deux lumières avec Blue FlexTip
- de** Zweilumiger, zentraler Venenkatheter mit Blue-FlexTip
- it** Catetere per vena centrale a due lumi con Blue FlexTip
- pt** Cateter Venoso Central com Lúmen Duplo e Blue FlexTip
- es** Catéter para vena central de dos luces con Blue FlexTip
- sv** Tvåkanalig centralvenös kateter med Blue FlexTip

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65
Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.57	5000
Proximal (18 Ga.)	0.40	1500

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IC-17702-135A (03/19)



Two-Lumen Central Venous Catheter with Blue FlexTip®

REF **MC-17702**

7 Fr. | **2** Lumen | **20** cm catheter length