



Two-Lumen Hemodialysis Catheter with Blue FlexTip® Arrowg+ard Blue® for High Volume Infusions

REF **MC-25122-F** 12 Fr. | 2 Lumen | 20 cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

[†]Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


Rx only.

Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels. The Arrowg+ard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

- fr** Cathéter pour hémodialyse à deux lumières avec Blue FlexTip Arrowg+ard-Blue pour les perfusions de grand volume
- de** Zweilumiger Hämodialysekatheter mit Blue-FlexTip-Arrowg+ard-Blue für Infusionen großer Flüssigkeitsmengen
- it** Catetere per emodialisi a due lumi con Blue FlexTip Arrowg+ard-Blue per volumi di infusione elevati
- pt** Cateter para Hemodiálise com Lúmen Duplo e Blue FlexTip, Arrowg+ard Blue, para Infusões de Grandes Volumes
- es** Catéter de dos luces para hemodiálisis con Blue FlexTip Arrowg+ard Blue para infusiones de gran volumen
- sv** Tvåkanalig hemodialyskatheter med Blue FlexTip Arrowg+ard Blue för infusioner med hög volym



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.60	200	100
		340	200
Proximal (12 Ga.)	1.50	200	100
		340	200

* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IC-25122-130D (10/19)



Two-Lumen Hemodialysis Catheter with Blue FlexTip® Arrowg+ard Blue® for High Volume Infusions

REF **MC-25122-F** 12 Fr. | 2 Lumen | 20 cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

[†]Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


Rx only.

Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels. The Arrowg+ard Blue antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

- fr** Cathéter pour hémodialyse à deux lumières avec Blue FlexTip Arrowg+ard-Blue pour les perfusions de grand volume
- de** Zweilumiger Hämodialysekatheter mit Blue-FlexTip-Arrowg+ard-Blue für Infusionen großer Flüssigkeitsmengen
- it** Catetere per emodialisi a due lumi con Blue FlexTip Arrowg+ard-Blue per volumi di infusione elevati
- pt** Cateter para Hemodiálise com Lúmen Duplo e Blue FlexTip, Arrowg+ard Blue, para Infusões de Grandes Volumes
- es** Catéter de dos luces para hemodiálisis con Blue FlexTip Arrowg+ard Blue para infusiones de gran volumen
- sv** Tvåkanalig hemodialyskatheter med Blue FlexTip Arrowg+ard Blue för infusioner med hög volym



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.60	200	100
		340	200
Proximal (12 Ga.)	1.50	200	100
		340	200

* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



CE REP Teleflex Medical
IDA Business and Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IC-25122-130D (10/19)