

REF REF BD-09880



Percutaneous Sheath Introducer Kit
with Integral Hemostasis Valve/Side Port for use with
7 - 7.5 Fr. Catheter

8.5^{Fr.} | 10^{cm} sheath length | .035^{inch dia.} spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 2: Syringes: 5 mL Luer-Slip
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Drape: 24" x 36" with 4" fenestration
- 1: CSR Wrap
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Obturator : 8 Fr.
- 1: Scalpel: #11
- 2: Pre Sponge Swabs
- 1: Suture: 3-0 Silk with Straight Needle
- 1: Tray: Prep
- 1: SharpsAway® Disposal Cup

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**
Sterilized using ethylene oxide unless otherwise indicated in the contents list.
Fluid path components are non-pyrogenic.
This product contains dry natural rubber.



STERILE EO



Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA
Teleflex®
www.teleflex.com
ID-09880-102D (2021-04)

REF BD-09880



Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use
with 7 - 7.5 Fr. Catheter

8.5^{Fr.} | 10^{cm} sheath length | .035^{inch dia.} spring-wire guide