

# Peripherally Inserted Central Catheter with Seldinger Access Set with 80 cm Platinum Tipped, Marked Wire



## Arrow Select Kits

REF PR-01351-MH

### Contents:

- 1: Indwelling Catheter, Trimmable: 3 Fr. x 21-1/2" (55 cm) Radiopaque Polyurethane with Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard
- 1: Sheath: 17 Ga. x 2-3/4" (7 cm) Radiopaque Peel-Away over 3 Fr. Dilator
- 1: Spring-Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 31-1/2" (80 cm) (Straight Soft Tip on One End – Straight Stiff Tip on Other)
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 2-3/4" (7 cm) TW with Echogenic Tip
- 1: Syringe: 5 mL Luer-Slip
- 1: Catheter Trimmer
- 1: CSR Wrap
- 1: Chart Sticker
- 1: Patient ID Card
- 1: Patient Information Booklet
- 1: Paper Tape Measure
- 1: Injection Site Cap

**Catheter Diameter:** 3 Fr., 0.99 mm

**Priming Volume: 0.30 mL**

Priming volumes are approximate and are done without the injection cap  
Injection cap priming volume is 0.17 mL.

**Flow Rate: 150 mL/hr**

Flow rates are done with normal saline, room temperature, 100 cm head height  
and represent approximate flow capabilities.

**Pump Flow Rate: 600 mL/hr**

Pump flow rates are determined at maximum pump pressure of 10 psig.



**Rx only.**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide.  
Do not use if package has been previously opened or damaged.  
Fluid path components are non-pyrogenic. Contains no medication.

Arrow International, Inc.  
2400 Bernville Road,  
Reading, PA 19605 USA

EC REP Teleflex Medical  
IDA Business and Technology Park,  
Athlone, Ireland



IH-01351-100B (5/16)

# Peripherally Inserted Central Catheter with Seldinger Access Set with 80 cm Platinum Tipped, Marked Wire

REF PR-01351-MH

Arrow Select Kits



3 Fr.

55 cm  
catheter  
length

.016/.018 inch diameter  
spring-wire  
guides