

Product Includes
One-Piece
Sheath
Design

Contents:


- 1: Sheath: 6 Fr. x 4" (10 cm) Arrow-Flex®
Radiopaque Polyurethane with Integral
Hemostasis Valve/Side Port and Tissue Dilator
1: Spring-Wire Guide: .035" (0.89 mm) dia.
x 17-3/4" (45 cm) (Straight Soft Tip
on One End – "J" Tip on Other)
1: Introducer Needle: 18 Ga. x 2-1/2"
(6.35 cm) XTW
1: Cath-Gard® Catheter Contamination Shield:
30.5 cm with Sterile Tape
1: Obturator: 8 Fr.
3: Gauze Pads: 4" x 4"

Rx only.


Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.

This Product Contains Dry Natural Rubber.

 Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

 Teleflex Medical
IDA Business and Technology Park
Athlone, Ireland



REF SI-09600

Percutaneous Sheath
Introducer Set
with Integral Hemostasis Valve/Side Port
for use with 4 - 6 Fr. Catheters

ARROW®

- FR** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 4 - 6 Fr.
- DE** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 4 - 6 Fr.
- IT** Set per introduzione di guaina percutanea con valvola emostatica/ bocchetta laterale incorporate per l'uso con cateteri da 4 - 6 Fr.
- PL** Zestaw do przezskórnego wprowadzania koszulki, z integralnym zaworem hemostatycznym / portem bocznym, do użycia z cewnikami o rozmiarze 4 - 6 F
- PT** Conjunto de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 4 - 6 Fr.
- RU** Комплект чрескожной капсулы-интубатора с несъемным гемостазным клапаном/боковым портом, для использования с катетерами 4 - 6 Fr.
- SL** Uvajalni komplet s perkutanim tulcem z vgrajenim hemostatskim ventilom/stranskim vhodom za uporabo s katetri 4 - 6 Fr.
- ES** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 4 a 6 Fr.
- SV** Set för perkutan skyddsintroducerare med integral hemostasventil/ sidoöppning för bruk med katetrar mellan 4 och 6 Fr.
- TR** Perkütan Kılıf Introducer Seti, Entegre Hemostaz Valfi/Yan Portlu, 4 - 6 Fr. Kateterler ile kullanılmak üzere

1109600-1308 (8/11)

ARROW®
INTERNATIONAL

REF
SI-09600

6^{Fr.}

10^{cm}
sheath
length

2.9^{mm}
sheath
O.D.

.035^{inch dia.}
spring-wire
guide

Percutaneous Sheath Introducer Set

Peel to
Open