



REF SI-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide

MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

Contents:

- 1: Two-Lumen Central Venous Access Device¹: 10 cm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 20.7 cm
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Tuohy-Borst Adapter, Sterile Tape
- 1: Dust Cap: Vented
- 3: Dust Cap: Non-Vented
- 1: Stopcock: 4-way High-Flow
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

¹Licensed under US Patent No. 5,250,038.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.


⚠ California Prop. 65

Warning: Cancer and Reproductive Harm
www.P65Warnings.ca.gov.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu d'accès veineux central à deux lumières MAC avec valve hémostatique intégrée, à utiliser avec les cathéters de 7 - 8 Fr.
- de** Einführbesteck für zweilumigen, MAC-zentralen Venenkatheter mit Hämostaseventil mit Kathetern der Größen 7 - 8 Fr.
- it** Set di accesso venoso centrale a due lumi MAC con valvola emostatica incorporata per l'uso con cateteri da 7 - 8 Fr.
- pl** Dwukanałowy zestaw MAC do uzyskiwania dostępu do żył centralnych, z integralnym zaworem hemostatycznym, do użycia z cewnikami o rozmiarze 7 - 8 F
- pt** Conjunto para Acesso Venoso Central MAC com Lúmen Duplo, Válvula de Hemostase Integral, para utilizar com Cateters de 7 - 8 Fr.
- ru** Комплект для катетеризации главных вен MAC с двухпросветным катетером и несъемным гемостазным клапаном, для использования с катетерами 7 - 8 Fr.
- sl** Komplet za centralni venski dostop MAC z dvema svetlinama z vgrajenim hemostatskim ventilom za uporabo z katetri 7 – 8 Fr.
- es** Conjunto para Acceso Venoso Central MAC con Lúmen Duplo, Válvula de Hemostase Integral, para utilizar con Cateters de 7 - 8 Fr.
- sv** MAC tvåkanaligt set för centralt ventillträde med integrerad hemostasventil för användning med 7 - 8 Fr katetrar
- tr** MAC İki Lümenli Santral Venöz Erişim Seti, Entegre Hemostaz Valfi, 7 - 8 Fr. Kateterler ile kullanılmak üzere



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
0086 Co. Westmeath, Ireland

Teleflex®
www.teleflex.com
II-11142-131A (01/19)



REF SI-11142

9 Fr. distal lumen | 2 Lumen | 10 cm device length | .035 inch dia. spring-wire guide

MAC™ Two-Lumen Central Venous Access Set with Integral Hemostasis Valve for use with 7 - 8 Fr. Catheters

